cBioPortal Tutorial #4: Virtual Studies

Create and save your own custom study comprised of samples from existing studies

Last update: December 23, 2021

Tutorial Objectives

- Explain what virtual studies are and how to create them
- Delineate the different features available when a user is or is not logged in
- Show how to log in

Virtual Study Overview

What is a Virtual Study?

A virtual study is a custom study comprised of samples from one or more existing studies.


Why is a virtual study useful?

The virtual study feature allows you to define a custom cohort of samples that fit your specific genomic or clinical criteria of interest. These samples can be a subset of the data available in an existing study, or result from the combination of multiple existing studies. This cohort of samples can then be queried or explored just like a traditional study, and can be returned to at a later date or shared with a collaborator.

In this tutorial, we will create a virtual study comprised of a subset of samples from two existing TCGA glioma studies.


We begin by selecting the studies of interest and viewing them in the Study Summary page.

Creating a virtual study: Selecting studies



Creating a virtual study: Study summary page


This is the same Study Summary page that you see for a single study, but now with data from 2 studies.



Creating a virtual study: Study summary page

Link to this page


In this example, we will additionally filter to samples with mutations in IDH1.




We are now ready to create our virtual study.

What you see will differ slightly depending on if you are logged in or not. The next two slides show what the virtual study button looks like in each scenario.

Virtual study button (not logged in)



Virtual study button (logged in)




on to view this new virtual study.

Whether or not you are logged in, let's continue

Virtual study: Study summary page

We can now see the study summary page for our newly created virtual study. You can now explore these samples or run a query, just like with a traditional study. Remember that while making this study we filtered down to IDH1 mutant samples, but the non-mutant samples were still present. Now, those samples without mutations are completely gone.




Logging in to cBioPortal:

If you use a private instance of cBioPortal, you might already go through an authentication process (through your university, for example). If so, there is nothing else you need to do.

But if you don't normally log in, then let's go through how this works.

Logging in



Logging in


We use Google and Microsoft for authentication. Sign in with any Google or Microsoft account.



Welcome to cBioPortal - sign in with your Google or Microsoft account to store your virtual studies and groups. This will allow you to access your studies and groups from any computer, and cBioPortal will also remember your study view charts preferences for each study (i.e. order of the charts, type of charts and visibility). Login is optional and not required to access any of the other features of cBioPortal.



Logging in



Virtual Study Authentication Questions

Is it necessary to log in to use virtual studies?

No. A user that has not logged in can create virtual studies and run queries in those studies (by using the query box on the study summary page). Links to virtual studies are permanent, so you can save the link on your computer and come back to it anytime, or share it with others.

If I do log in, what additional functionality do I gain?

If you log in, you gain the ability to save your virtual study to the list of existing studies on the homepage. This makes a virtual study functionally the same as any other study: you can access your virtual studies in the query builder and you can combine an existing virtual study with any other study to create a new virtual study.

or email us at:

Questions?

Check out our other tutorials

cbioportal@googlegroups.com